

NEW

EndoVue® Plus 19"

Value Endoscopy Visualization



- Optimized for endoscopy applications
- Exceptional image quality
- Multi-modality imaging capabilities
- Low power consumption and fanless cooling
- HD-SDI and DVI inputs
- Made in the USA

For the medical professional seeking to maximize their capital investment and upgrade their aging monitor, the EndoVue® Plus series delivers the image performance of a high-end endoscopy display at a value price.

Manufactured in the USA and offering a 3-year warranty, EndoVue Plus is a versatile visualization display, supporting multi-modality medical imaging requirements. It is plug-and-play compatible with endoscopic camera systems from leading manufacturers. Designed and built for medical use, it is FDA listed, TUV certified, and CE marked, offering patient safety, clinical performance, and durability compared to non-medical grade displays.

The EndoVue Plus 19" display accommodates standard analog and digital HD-SDI and DVI-I high-definition signals from a variety of medical imaging sources, including ultrasound, PACS, and patient monitor vital signs. It features a fanless cooling design to minimize the risk of spreading airborne contaminants within the sterile field, and industry-leading LED backlight technology for reduced power consumption.

EndoVue® Plus 19"

Value Endoscopy Visualization



SPECIFICATIONS

| | |
|-----------------------|--|
| Screen Size | 19" diagonal |
| Image Size | 14.8" x 11.9"/376 x 301 mm (WxH) |
| Resolution | 1280 x 1024 |
| Aspect Ratio | 5:4 |
| View Angle | 170° |
| Luminance | 350 cd/m ² typical |
| Contrast Ratio | 1500:1 (typical) |
| Numbers of Colors | 16.8 Million (8-bit) |
| Response Time | 15ms |
| Latency | Up to 16ms |
| Video Inputs | HD-SDI/DVI-I/RGHS/VGA/YPbPr/S-Video/Composite |
| Serial Control | RS-232 (DB-9) |
| DC Output | 5V (1.0 Amp max) |
| Warranty | 3 years |
| Dimensions | 18" x 16" x 4"/465 x 400 x 98 mm (WxHxD) |
| Weight | 9.0 lbs/4.2 kg |
| Shipping Weight* | 21 lbs/9.5 kg |
| Shipping Dimensions | 24" x 8" x 23"/61 x 20 x 58 cm |
| Green Compliance | REACH, RoHS-2, WEEE, Conflict Minerals |
| Regulatory Compliance | ANSI/AAMI ES60601-1, CAN/CSA C22.2 No. 60601-1, CAN ICES-3(B)/NMB-3(B), FCC Class B, EN60601-1, EN60601-1-2, CE, MDD 93/42/EEC, CCC, Class I Medical Device, EAC |
| Supplied Accessories | Power Supply, 100-240 VAC, 50-60 Hz to 24 VDC, 2.7 A AC power cord, region appropriate User Manual – English with multi-language CD-ROM |

*Shipping weights are approximate.

ORDER INFORMATION

EndoVue Plus 19

| | |
|---------|------------------|
| 90K0016 | EndoVue Plus 19" |
|---------|------------------|

OPTIONAL ACCESSORIES

| | |
|---------|-----------------------------------|
| 26B0102 | Display Stand - 19", 21", 24" |
| 35X0100 | Fiber Converter 5 VDC Power Cable |

Global Headquarters

www.ndssi.com
5750 Hellyer Ave
San Jose, CA 95138
USA 🇺🇸
+1 408 776 0085
info@ndssi.com

Europe

Novanta Europe GmbH
Münchner Strasse 2A
82152 Planegg
Germany 
+49 89 31 707 100
info.novanta-europe@novanta.com

Asia Pacific

Novanta Japan
East Square Omori
6-20-14 Minamioi, Shinagawa-ku
Tokyo 140-0013
Japan
+81 3 5753 2466
info@ndssi.jp

Oriental Media Center
Suite 2302, Tower C
No.4, Guang Hua Road
Chao Yang District
Beijing, 100026
China
+86 10 8559 7859
info@ndssi.com.cn

NDS Quality System

ISO 9001:2008 and ISO 13485:2003
FDA Registration #2954921