

September 14, 2016

To Our Valued Customers:

NDS announces the **End-of-Life (EOL) and Last Time Buy (LTB)** of the **ZeroWire Ultra** effective **December 31, 2016 or while supplies last**. After this date, this product will no longer be available for purchase and **current inventory is limited**.

NDS encourages customers to transition to the suggested replacement products listed below. To fulfill our customers' needs, there will be last time buy (LTB) opportunity. Please submit your last PO by September 30, 2016. Supplies are limited, and are available on a first-come, first-serve basis.

End-of-Life Product

- 90T2035 ZeroWire Ultra - Duo, DVI
- 90T2036 ZeroWire Ultra - Rx, DVI
- 90T2037 ZeroWire Ultra - Tx, DVI
- 90T2038 ZeroWire Ultra - Duo, DVI, SDI
- 90T2039 ZeroWire Ultra - Tx, DVI, SDI
- 90T2040 ZeroWire Ultra - Duo, DVI, Single Channel
- 90T2041 ZeroWire Ultra - Rx, DVI, Single Channel
- 90T2042 ZeroWire Ultra - Tx, DVI, Single Channel
- 90T2043 ZeroWire Ultra - Duo, DVI, SDI, Single Channel, China
- 90T2044 ZeroWire Ultra - Rx, DVI, Single Channel, China
- 90T2045 ZeroWire Ultra - Tx, DVI, SDI, Single Channel, China
- 90T2060 ZeroWire Ultra - Duo, DVI, SDI, Single Channel
- 90T2061 ZeroWire Ultra - Tx, DVI, SDI, Single Channel
- 90T2062 ZeroWire Ultra - Rx, DVI
- 90T2063 ZeroWire Ultra - Tx, DVI
- 90T2064 ZeroWire Ultra - Tx, DVI, SD
- 90T2004 ZeroWire Classic - Duo, DVI
- 90T2005 ZeroWire Classic - Tx, DVI
- 90T2006 ZeroWire Classic - Rx, DVI
- 90T2007 ZeroWire Classic - Duo, DVI, SDI
- 90T2008 ZeroWire Classic - Tx, DVI, SDI
- 90T2014 ZeroWire Classic - Rx, DVI, Accessories
- 90T2015 ZeroWire Classic - Tx, DVI, Accessories
- 90T2016 ZeroWire Classic - Tx, DVI, SDI, Accessories, Single Channel
- 90T2024 ZeroWire Classic - Duo, DVI, Single Channel
- 90T2032 ZeroWire Classic - Rx, DVI, Accessories, Single Channel, China
- 90T2033 ZeroWire Classic - Tx, DVI, SDI, Accessories, Single Channel, China
- 90T2034 ZeroWire Classic - Duo, DVI, SDI, Accessories, Single Channel, China
- 90T2058 ZeroWire Classic - Duo, DVI, SDI, Single Channel
- 90T2059 ZeroWire Classic - Tx, DVI, SDI, Accessories, Single Channel

The next generation ZeroWire G2 will be functionally equivalent while providing improved performance.

Successor Product

- 90T2070 ZeroWire G2 Pair
- 90T2071 ZeroWire G2 Receive (Rx) Only
- 90T2072 ZeroWire G2 Transmit (Tx) Only
- 90T2082 ZeroWire G2 Pair Japan version
- 90T2083 ZeroWire G2 Receive (Rx) Only Japan version
- 90T2084 ZeroWire G2 Transmit (Tx) Only Japan version

Please refer to the new ZeroWire G2 Sales Sheet and User Manual at www.ndssi.com or contact your local sales manager or regional office for details regarding accessory kits, specifications, features and benefits of this new HD video wireless system:

- Uncompressed 60 GHz technology provides superior image quality
- Small footprint, all-in-one wireless solution with built-in 1080p scaling
- 2D and 3D video support (refer to Sales Sheet for supported 3D formats)
- Real-time feedback via on-screen overlays for status and signal strength
- Advanced Auto-Tracking system with an angled antenna system

We would like to take this opportunity to express our great appreciation for choosing NDS as your display partner. Our customer's voice is valued and always welcome. We are committed to supporting your transition to the new ZeroWire G2. Should you have any questions or would like to receive any further information, please do not hesitate to let us know.

Sincerely,

NDS Sales and Marketing Teams

- Corporate/Americas: +408 776 0085 info@ndssi.com
- EMEA: +49 89 31 707 100 info.gsig-europe@gsig.com
- China: +86 10 8559 7859 info@ndssi.com.cn
- Japan: +81 3 5753 2466 info@ndssi.jp

