

LIFEVUE™ SERIES

Patient Monitoring Displays



LifeVue



The LifeVue™ series of displays from NDSsi are designed specifically for patient monitoring applications, capable of rendering vital sign wave forms without image-lag artifacts. These fully customizable display solutions are the highest quality and most versatile patient monitoring displays available today. Fully compliant for medical use, LifeVue™ displays offer exceptional image clarity, reliability, and safe operation. Available with optional touch screen and audio capabilities, these displays are compatible with patient monitoring equipment from all major manufacturers and are designed to meet the medical requirements for respiratory care signal alarms.

- Medical Grade 15" LED and 19" Sizes with Touch and Audio
- Industry Leading LED Backlight Technology*
- Fully Compatible with all Monitoring Equipment
- Easy-To-Use and Customize (Optional Interfaces and Features)
- OR/ER/ICU Patient Monitoring

First in Imaging

* Applies to 15" LifeVue product only.



Touch Screen



19"

15" LED

TYPICAL APPLICATIONS

- Bedside patient monitoring
- Centralized care station monitoring
- Integrated OR touch panel control interface
- OR/ER environments for patient monitoring
- Mobile cart mounted systems

PRODUCT DETAILS

- 15" LED or 19" diagonal size
- Touch screen and audio alarm support
- Advanced security: control interface lockout
- 100 mm VESA mounting
- Optional desktop stands available

LIQUID INGRESS

- Designed to meet liquid ingress standard EN-60529 (IEC-529)

INTERFACE

- Touch screen for a mouse-less environment
- Accepts both DVI-Digital or DVI-Analog (VGA) video inputs
- Supports standard VESA resolution graphics input signals

RESPIRATORY AUDIO ALARM

- Provides >70dB of sound at 1 meter in front of the display
- Designed to meet medical requirements of ISO/DIS 9703-2 Respiratory Care Signal Alarms

NDSSI QUALITY SYSTEM

ISO 9001:2008 and ISO 13485:2003; FDA Registration # 2954921

COMPLIANCE & CERTIFICATIONS

ANSI/AAMI ES60601-1, CAN/CSA C22.2 No. 60601-1, FCC Class B, EN60601-1, EN60601-1-2, CE, MDD 93/42/EEC, 2007/47/EC, RoHS, WEEE, Class I Medical Device, IP32 (Pending)



Corporate Headquarters

5750 Hellyer Avenue
San Jose, CA 95138 (USA)
Tel: 408 776 0085
Toll Free: 866 637 5237
Email: info@ndssi.com

Europe

Nijverheidscentrum 28
2761 JP Zevenhuizen (ZH)
The Netherlands
Tel: + 31 180 63 43 56
Email: info-EMEA@ndssi.com

Asia Pacific

East Square Omori
6-20-14 Minamioi, Shinagawa-ku
Tokyo 140-0013, Japan
Tel: +81 3 5753 2466
Email: info@ndssi.jp