

NEW

## EndoVue® 24"

### Value Endoscopy Visualization



## HD ENDOVIEW

- Optimized for endoscopy applications
- Outstanding image quality
- Multi-modality imaging capabilities
- Low power consumption and fanless cooling
- Robust and reliable
- Minimum cost of ownership

For the medical professional seeking to maximize their capital investment and upgrade their aging monitor, the EndoVue® series delivers the image performance of a high-end endoscopy display at a value price.

EndoVue is a versatile visualization platform, supporting multi-modality medical imaging requirements. It is plug-and-play compatible with endoscopic camera systems from leading manufacturers. Designed and built for medical use, it is FDA listed, TUV certified, and CE marked, offering patient safety, clinical performance, and durability compared to non-medical grade displays.

The EndoVue display accommodates standard analog and digital high-definition signals from a variety of medical imaging sources, including ultrasound, PACS, and patient monitor vital signs. It features a fanless cooling design to minimize the risk of spreading airborne contaminants within the sterile field, and industry-leading LED backlight technology for reduced power consumption.

# EndoVue® 24"

## Value Endoscopy Visualization

### SPECIFICATIONS

Screen Size	24" diagonal
Image Size	20.4" x 12.8"/518 x 324 mm (WxH)
Resolution	1920 x 1200
Aspect Ratio	16:10
View Angle	178°
Luminance	300 cd/m <sup>2</sup> typical
Contrast Ratio	1,000:1 (typical)
Numbers of Colors	16.7 Million (8-bit)
Response Time	14ms
Latency	20ms
Video Inputs	DVI-D, RGBS/VGA/YpPr (HD-15 or DVI-I via adapter), S-Video, Composite
Video Output	DVI-D
Warranty	1 Year standard, optional extended 1 year
Dimensions	21.9" x 16.1" x 2.8"/566 x 408 x 70 mm (WxHxD)
Weight	14.7 lbs/6.7 kg
Shipping Weight*	22 lbs/9.9 kg
Shipping Dimensions	25" x 8" x 23"/64 x 21 x 59 cm
Green Compliance	REACH, RoHS-2, WEEE, Conflict Minerals
Regulatory Compliance	ANSI/AAMI ES60601-1, CAN/CSA C22.2 No. 60601-1, CAN ICES-3(B)/NMB-3(B), FCC Class B, EN60601-1, EN60601-1-2, CE, MDD 93/42/EEC, CCC, EAC, Class I Medical Device
Supplied Accessories	Power Supply, 100-240 VAC, 50-60 Hz to 24 VDC, 2.7 A AC power cord, region appropriate User Manual – English with multi-language CD-ROM

\*Shipping weights are approximate.

### ORDER INFORMATION

#### EndoVue 24

90K0070	EndoVue 24" - North America
90K0071	EndoVue 24" - South America
90K0072	EndoVue 24" - Japan
90K0073	EndoVue 24" - China
90K0074	EndoVue 24" - Emerging Market**
90K0075	EndoVue 24" - Europe

### OPTIONAL ACCESSORIES

26B0102	Display Stand - 19", 21", 24"
35X0100	Fiber Converter 5 VDC Power Cable
73A1015	EndoVue 24" Extended 1 Year Warranty

\*\*Emerging Market model eligible countries (non-eligible countries will void the warranty)

Bahrain | Bangladesh | Botswana | Bulgaria | Côte d'Ivoire | Croatia | Cyprus | Czech Republic | Egypt | Estonia | Ghana | Hungary | India | Indonesia | Jordan | Kenya | Lithuania | Macedonia | Malaysia | Malta | Mauritius | Morocco | Nigeria | Oman | Pakistan | Philippines | Poland | Qatar | Romania | Russia | Serbia | Slovakia | Slovenia | South Africa | Sri Lanka Taiwan | Thailand | Tunisia | Turkey | UAE | Vietnam



#### Global Headquarters

www.ndssi.com  
5750 Hellyer Ave  
San Jose, CA 95138  
USA 🇺🇸  
+1 408 776 0085  
info@ndssi.com

#### Europe

Novanta Europe GmbH  
Münchner Strasse 2A  
82152 Planegg  
Germany   
+49 89 31 707 100  
info.novanta-europe@novanta.com

#### Asia Pacific

Novanta Japan  
East Square Omori  
6-20-14 Minamioi, Shinagawa-ku  
Tokyo 140-0013  
Japan  
+81 3 5753 2466  
info@ndssi.jp

Oriental Media Center  
Suite 2302, Tower C  
No.4, Guang Hua Road  
Chao Yang District  
Beijing, 100026  
China  
+86 10 8559 7859  
info@ndssi.com.cn

#### NDS Quality System

ISO 9001:2008 and ISO 13485:2003  
FDA Registration #2954921